

Ersatzlampen für den Dentalbereich
Replacement lamps for Dental use

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X446		„K“ Super Starlite	24	150	GY9.5
99-X449	78-8045-4571-9	3 M Ortholux / Visolux I / Visolux II	12	75	G5.3-4.8
99-X449	78-8045-4571-9	3 M XL 1000 / XL 3000	12	75	G5.3-4.8
99-X453		Acta	14	35	GZ4
99-X421		Algotec Walklite	8	50	GZ6.35
99-X424		Algotech LH 5000	15	150	GZ6.35
99-X449	3442	Almore OP II	12	75	G5.3-4.8
99-X454		Alpha 1975	6,6	200	R7s
99-X446		Amaca Elisse	24	150	GY9.5
99-X461		Arcus I / ABL-5	6	35	GZX4
99-X462		Arcus II	10	52	
99-X449		Astral	12	75	G5.3-4.8
99-X449		Bayer	12	75	G5.3-4.8
99-X450		Celebrity	12	55	PK22s
99-X450		Chayes Virginia	12	55	PK22s
99-X424	4070110	Coe-Lite	15	150	GZ6.35
99-X449	C7999	Coltolux II / 3 / 4	12	75	G5.3-4.8
99-X451		Cromalux 100	12	100	G5.3-4.8
99-X401	SYC0023	Cuda Products M3-150 Multi Light source	21	150	GX7.9
99-X455		Cygna-Scope	120	250	GY5.3
99-X437		D 44 Neu	24	150	G6.35
99-X449	950-9022	Darby Superdent 2001	12	75	G5.3-4.8
99-X456	950-9024	Darby Superdent 7500	21	150	GX5.3
99-X449		De Trey Prismetics	12	75	G5.3-4.8
99-X451		De Trey Prismetics Lite	12	100	G5.3-4.8
99-X451		De Trey Prismetics Lite 2	12	100	G5.3-4.8
99-X449		Degulux	12	75	G5.3-4.8
99-X449	20437	Demetron Optilux 101 / 400 / 401	12	75	G5.3-4.8
99-X451	20140	Demetron Optilux 400 High Output	12	100	G5.3-4.8
99-X451	20140	Demetron Optilux 401 High Output	12	100	G5.3-4.8
99-X449	20437	Demetron VCL 200 / 300	12	75	G5.3-4.8
99-X456	3904	Den-Mat Visar II	21	150	GX5.3
99-X457	70353	Dentsply / York Triad 2000	120	250	GY5.3
99-X451		Dentsply Prismetics Lite / Lite 2	12	100	G5.3-4.8
99-X449		Dexter Destilux	12	75	G5.3-4.8
99-X456		Dolan Jenner Fiberlite 170 / 176	21	150	GX5.3
99-X451	3021	Efos-Lite	12	100	G5.3-4.8
99-X449		Eggert-Lux	12	75	G5.3-4.8
99-X449		Elipar 2	12	75	G5.3-4.8
99-X426		Euro 90	13,8	50	GX5.3
99-X437		Eurodent	24	150	G6.35
99-X430		Eurodent	8	20	GZX4
99-X453		Euromax	14	35	GZ4
99-X454		Fantastic	6,6	200	R7s
99-X437		Faro S70	24	150	G6.35
99-X430		Flex	8	20	GZX4
99-X446		G208 (1. Version)	24	150	GY9.5
99-X450		Healthco Aristocrate	12	55	PK22s
99-X449	HCOB400491	Healthco Aristocrate I	12	75	G5.3-4.8
99-X453		Healthco Aristocrate II	14	35	GZ4
99-X424	HCOB400711	Healthco High Intensity VLC Unit	15	150	GZ6.35
99-X457	HCOB400480	Healthco VLC Illuminator (Wand)	120	250	GY5.3
99-X423		Heliomat	12	100	GZ6.35
99-X451		Heliospad	12	100	G5.3-4.8
99-X451	549-3519	Henry Schein VCL 300	12	100	G5.3-4.8
99-X437		Isolight	12	150	G6.35

Ersatzlampen für den Dentalbereich
Replacement lamps for Dental use

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X422		Johnson & Johnson Fotofill	12	75	GZ6.35
99-X453		Kavo H2	14	35	GZ4
99-X457	555-5001	Kavo Vicon DLS	120	250	GY5.3
99-X457	555-5001	Kavo Vicon LS Illuminator	120	250	GY5.3
99-X457	555-5001	Kavo Vicon Polar II	120	250	GY5.3
99-X463	5533881	Kavo Multiflex Kupplungen und Motoren	3,3	715 mA	-
99-X437		Kavolux (alt)	24	150	G6.35
99-X441		Kavolux (neu)	24	150	G6.35
99-X453		Kerr 50 / 150 / 250 / 350	14	35	GZ4
99-X449		Kerr 101 / 400 / 401 / 402	12	75	G5.3-4.8
99-X451		Kerr 400 / 401 / 402	12	100	G5.3-4.8
99-X458	15512	Kerr Manufacturing Command Light	21	150	GX5.3
99-X449		Kerr VCL200 / VCL 300	12	75	G5.3-4.8
99-X456	1450	Kinetic Instruments Sunlite II	21	150	GX5.3
99-X451	1511	Kinetic Instruments Sunspot	12	100	G5.3-4.8
99-X424		Kulzer Translux	15	150	GZ6.35
99-X451	2244	Kulzer Translux CL	12	100	G5.3-4.8
99-X424	2045	Kulzer Translux EC	15	150	GZ6.35
99-X412	644101	L.D. Caulk PR 1 / Prisma Lite	24	200	GX5.3
99-X424		L.D. Caulk Prisma Lite	15	150	GZ6.35
99-X451	644358	L.D. Caulk Prismetics Lite	12	100	G5.3-4.8
99-X449		L.D. Caulk Prismetics Lite II	12	75	G5.3-4.8
99-X457	10037	Lares Quartz Fire (SX)	120	250	GY5.3
99-X430		LC	8	20	GZX4
99-X446		Lineo o Scorpio	24	150	GY9.5
99-X424		Litema	15	150	GZ6.35
99-X449		Litex 660	12	75	G5.3-4.8
99-X459		Luxomat	230	9	G23
99-X430		Lysta	8	20	GZX4
99-X449		M 1	12	75	G5.3-4.8
99-X453		M 1	14	35	GZ4
99-X449		Megalux CS	12	75	G5.3-4.8
99-X438		Mega-Physik Dentalux Optima EC250	24	250	G6.35
99-X451		Merz HL2000	12	100	G5.3-4.8
99-X446		Midway	12	150	GY9.5
99-X457	490874	Midwest Insight I / II	120	250	GY9.5
99-X460	490709	Midwest Oral Illuminator II / III	21	150	GX7.9
99-X444		Mikrofilmlesegerät	13.8	85	GX5.3
99-X450		Moonlight	12	55	PK22s
99-X423		NSK Opticure	12	100	GZ6.35
99-X437		OMS	24	150	G6.35
99-X440		OP Siemens	12	100	GY6.35
99-X449		Optilux UV	12	75	G5.3-4.8
99-X402		O'Ryan Industries Omega 110 / 510	120	150	GX7.9
99-X456	301203	Penn-Med 200L (Wand)	21	150	GX5.3
99-X449	301315	Penn-Med 310	12	75	G5.3-4.8
99-X449	3050	Phasealloy Whitelite	12	75	G5.3-4.8
99-X424		Pluraflex HL 150	15	150	GZ6.35
99-X453		Polilux Kaltl.	14	35	GZ4
99-X449		Polofil / Pololux	12	75	G5.3-4.8
99-X459		Polygerät	230	9	G23
99-X430		Polylight	8	20	GZX4
99-X421		Polysteril / Polysteril II	8	50	GZ6.35
99-X424		PR Lux 200	15	150	GZ6.35
99-X451	20-055	Pro Den Systems Executor	12	100	G5.3-4.8
99-X451		Profilux	12	100	G5.3-4.8

Ersatzlampen für den Dentalbereich

Replacement lamps for Dental use

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X453		ProLite	14	35	GZ4
99-X451	17244	S.S. White Profile VCS	10	100	G5.3-4.8
99-X446		S2000	24	150	GY9.5
99-X421		Satelec Servolite	8	50	GZ6.35
99-X451	64624	Schutz Dental Elcelux	12	100	G5.3-4.8
99-X457		Sci Can Dual-Lux II (Curing)	120	250	GY5.3
99-X452		Sci Can Dual-Lux II / III(Illum.)	21	150	GX5.3
99-X411		Sci Can Dual-Lux III (Curing)	21	150	GX5.3
99-X424		Siemens C1	15	150	GZ6.35
99-X464	5940291	Siemens, Lampe für Kupplungen	-	-	-
99-X554	5945407	Siemens, Lampe für Motoren	-	-	-
99-X450		Sirolux	12	55	PK22s
99-X440		Sirolux S	12	100	GY6.35
99-X437		Sirolux S	24	150	G6.35
99-X421		SL 4055	8	50	GZ6.35
99-X424		Solid State Systems Initiator	12	150	GZ6.35
99-X401	61351	Star Dental Starlite	21	150	GX7.9
99-X446		Star Track	24	150	GY9.5
99-X437		Starlite (neu)	24	150	G6.35
99-X446		Stern Multilux	24	150	GY9.5
99-X431		Streppel Leuchte	6	20	G4
99-X446		Sunlight S90	24	150	GY9.5
99-X451	41501	Teledyne Getz Focus Activator	12	100	G5.3-4.8
99-X437		Ultradent (alt)	24	150	G6.35
99-X449		Visilux / Vilolux II	12	75	G5.3-4.8
99-X449		Vivadent DLX	12	75	G5.3-4.8
99-X411	7210	Vivadent Heliomat	21	150	GX5.3
99-X424	7110	Vivadent Heliomat	155	150	GZ6.35
99-X453		Vivadent Heliolux / Heliolux II	14	35	GZ4
99-X449		Vivadent Heliolux GTE	12	75	G5.3-4.8
99-X461		Vivalux	6	35	GZX4
99-X449		Vivalux II	12	75	G5.3-4.8

Metalldampflampen mit Reflektor, u.a. für Lichtquellen

Metal halide lamps with mirror

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X756	HCI 400W		80 V	400 W	
99-X752	HMI 575W/GS		100 V	575 W	
99-X740	HMP 575 SE		100 V	575 W	G 22
99-X742	MSR 575 HR		95 V	575 W	G 22
99-X736	HSR-400		67 V	400 W	GX9,5
99-X232SE	HTI 250 W/SE	z.B. Storz 20123120	45 V	270 W	FaX 1.5
99-X222	HTI 250 W/22		45 V	270 W	Spezial
99-X222C	HTI 250 W/22 C		45 V	270 W	Litzen + Stecker
99-X232	HTI 250 W/32		45 V	270 W	Spezial
99-X232C	HTI 250/32 C		45 V	270 W	Litzen + Stecker
99-X200	HTI 400 W/24		55 V	400 W	Spezial
99-X755	HTI 403 W/24		55 V	400 W	Litzen + Stecker
99-X744	HQI-TS150/NDL			150 W	R7s
99-X739	VIP R 273/45		38 V	270 W	Spezial
99-X751	VIP R 400/32			400 W	Spezial
99-X335	GE Gemini 300	Ansi EZG	35 V	300 W	Stecker
99-X3375	GE Marc 300/16	Ansi EZM	37.5 V	300 W	Stecker
99-X345	GE Marc 350/16T	Ansi EZT	45V	350W	Stecker