

Ersatzlampen für den Dentalbereich

Replacement lamps for Dental use

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X453		ProLite	14	35	GZ4
99-X451	17244	S.S. White Profile VCS	10	100	G5.3-4.8
99-X446		S2000	24	150	GY9.5
99-X421		Satelec Servolite	8	50	GZ6.35
99-X451	64624	Schutz Dental Elcelux	12	100	G5.3-4.8
99-X457		Sci Can Dual-Lux II (Curing)	120	250	GY5.3
99-X452		Sci Can Dual-Lux II / III(Illum.)	21	150	GX5.3
99-X411		Sci Can Dual-Lux III (Curing)	21	150	GX5.3
99-X424		Siemens C1	15	150	GZ6.35
99-X464	5940291	Siemens, Lampe für Kupplungen	-	-	-
99-X554	5945407	Siemens, Lampe für Motoren	-	-	-
99-X450		Sirolux	12	55	PK22s
99-X440		Sirolux S	12	100	GY6.35
99-X437		Sirolux S	24	150	G6.35
99-X421		SL 4055	8	50	GZ6.35
99-X424		Solid State Systems Initiator	12	150	GZ6.35
99-X401	61351	Star Dental Starlite	21	150	GX7.9
99-X446		Star Track	24	150	GY9.5
99-X437		Starlite (neu)	24	150	G6.35
99-X446		Stern Multilux	24	150	GY9.5
99-X431		Streppel Leuchte	6	20	G4
99-X446		Sunlight S90	24	150	GY9.5
99-X451	41501	Teledyne Getz Focus Activator	12	100	G5.3-4.8
99-X437		Ultradent (alt)	24	150	G6.35
99-X449		Visilux / Vilolux II	12	75	G5.3-4.8
99-X449		Vivadent DLX	12	75	G5.3-4.8
99-X411	7210	Vivadent Heliomat	21	150	GX5.3
99-X424	7110	Vivadent Heliomat	155	150	GZ6.35
99-X453		Vivadent Heliolux / Heliolux II	14	35	GZ4
99-X449		Vivadent Heliolux GTE	12	75	G5.3-4.8
99-X461		Vivalux	6	35	GZX4
99-X449		Vivalux II	12	75	G5.3-4.8

Metalldampflampen mit Reflektor, u.a. für Lichtquellen

Metal halide lamps with mirror

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-X756	HCI 400W		80 V	400 W	
99-X752	HMI 575W/GS		100 V	575 W	
99-X740	HMP 575 SE		100 V	575 W	G 22
99-X742	MSR 575 HR		95 V	575 W	G 22
99-X736	HSR-400		67 V	400 W	GX9,5
99-X232SE	HTI 250 W/SE	z.B. Storz 20123120	45 V	270 W	FaX 1.5
99-X222	HTI 250 W/22		45 V	270 W	Spezial
99-X222C	HTI 250 W/22 C		45 V	270 W	Litzen + Stecker
99-X232	HTI 250 W/32		45 V	270 W	Spezial
99-X232C	HTI 250/32 C		45 V	270 W	Litzen + Stecker
99-X200	HTI 400 W/24		55 V	400 W	Spezial
99-X755	HTI 403 W/24		55 V	400 W	Litzen + Stecker
99-X744	HQI-TS150/NDL			150 W	R7s
99-X739	VIP R 273/45		38 V	270 W	Spezial
99-X751	VIP R 400/32			400 W	Spezial
99-X335	GE Gemini 300	Ansi EZG	35 V	300 W	Stecker
99-X3375	GE Marc 300/16	Ansi EZM	37.5 V	300 W	Stecker
99-X345	GE Marc 350/16T	Ansi EZT	45V	350W	Stecker