


Laryngoskope

Laryngoscopes



Baby- und Kleinkindlaryngoskop  
Laryngoskop-Sets

Baby Laryngoscopes  
Sets of Laryngoscopes



## Baby-Laryngoskop

Bei der Entwicklung dieses handlichen Kleinstkind-Laryngoskopes wurden die anatomischen Verhältnisse des Neugeborenen besonders berücksichtigt. Störende aufgelötete oder abnehmbare Lampenträger sind hier nicht vorhanden. Die mit wenigen Handgriffen auswechselbaren Spatelblätter sind mit Längsrillen versehen. Diese ermöglichen eine gute Längsführung und verhindern ein seitliches Abrutschen des Spatels von der Zunge. Voll autoklavierbar (ohne Batterien) und dauerhaft resistent gegen Schaltstörungen durch bewussten Verzicht auf die Verwendung aufwendiger Schaltechnik. Ein- bzw. Ausschalten durch Herein- bzw. Herausdrehen der Bodenkappe.

## Baby-Laryngoscope

A handy set, especially constructed for new-borns, with regard for their anatomic conditions. Easy to clean, autoclaveable (without batteries). No disfunction because of broken switch or soldered connections; for switching on or off just screw the bottom cap in or out.

## Saling Laryngoskop-Set mit 3 Spateln

### 03-111

Batteriegriff (2x Typ AA), ø 19 mm, schmaler Spatel, breiter Spatel, Mikrospatel, 1 Glühlampe

## Saling-Laryngoscope Set with 3 blades

Battery handle (2x type AA), ø 19 mm, blade narrow, blade broad micro blade, 1 bulb



03-111

## Saling Laryngoskop-Set mit 2 Spateln

### 03-101

Batteriegriff (2x Typ AA), ø 19 mm, schmaler Spatel, breiter Spatel, 1 Glühlampe

## Saling-Laryngoscope Set with 2 blades

Battery handle, ø 19 mm, blade narrow, blade broad, 1 bulb

## Spatel für Saling-Laryngoskop

- ▶ **03-911** Spatel schmal
- ▶ **03-912** Spatel breit
- ▶ **03-913** Mikrospatel

## Blades for Saling-Laryngoscope

- ▶ Blade narrow
- ▶ Blade broad
- ▶ Micro blade



03-911



03-912



03-913

**Einzelteile für  
Saling-Laryngoskop**
**Components for  
Saling-Laryngoscope**

▶ <b>03-903</b> Batteriegriff (2x Typ AA), ø 19 mm	Battery handle (2x type AA), ø 19 mm,
▶ <b>99-006</b> Ersatzlampe, 2,5 V	Spare bulb, 2.5 V

