



# High power

LED-Stirnleuchten

LED headlights



www.faromed.de



## Basic-line .....

## Fantastic light

#### **▶** 08-795-S

LED head light, high precision light spot, wireless, basic model with **built in battery**, including charger 10-10795



White headband





#### **10-10795**

Charger



## Easy and quickly battery change







Easy and fast battery change

#### **▶** 08-797-S

LED head light, high precision light spot, wireless, basic model with easily **interchangeable battery**, including charger 10-10795



LiPo battery pack, 1x included



#### **10-10815**

Charger for 2x LiPo battery Option; not included





## Comfort-line.....

Touch dimmer, switch off automatic, robust headband for short time use (e.g. ENT)



#### Automatic On/Off Switch







Touch dimmer Easy hygenic cleaning, no knob



Easy and fast battery change

#### ▶ 08-796-S-1

LED head light, best for short time use (e.g. ENT), high precision light spot, wireless, basic model with 1x changeable battery, including charger 10-10795

## Smooth Headband

#### ► 08-806-S-1

Same as 08-796-S-1 but smooth headband for long procedures, touch dimmer, switch off automatic



#### **▶** 08-205-815

LiPo battery pack, 1x included

#### **▶ 10-10815**

Charger for 2x LiPo battery Option; not included



## High-end-line .....

## Extreme Powerful Light

#### **▶** 08-815

High performance, focusable headlight with High-Power-LED. Wireless, light weight, comfortable headlight with two interchangeable LiPo battery packs at the headband for working up to 7 hours and more.

Adjustable light intensity by "Touch dimmer"



#### **▶ 10-10815**

Charger for 2x LiPo battery



#### **▶** 08-205-815

LiPo battery pack, 2x included



## Light weight and high performance

#### **▶** 08-825

High performance, focusable headlight with High-Power-LED; no batteries at the headband for lowest weight and highest comfort. LiPo battery pack with belt clip, adjustable light intensity by "Touch dimmer"



#### **10-10780**

Charger for LiPo battery pack

#### **▶** 08-205-825

LiPo battery pack, 1x included

