

Xenon Headlights & Supplies

2011



► **Item code: 08-595-S**

Lightweight fiber optic headlight with high performance optic system. Comfortable and foam bolstered headlight (grey and black), easy adjustable with two central adjustment wheels. **Super lightweight cable**, l=2.5m, with gown clips, especially for long surgeries.



► **Item code: 08-590-S**

Lightweight fiber optic headlight with high performance optic system. Comfortable and foam bolstered headlight (grey and black); easy adjustable with central adjustment wheel. **Super lightweight cable**, l=2.5m, with gown clips, especially for long surgeries.



► **Item code: 08-585-S**

Comfortable fiber optic headlight with high performance optic system. **Sturdy, long-life fiber optic cable**, l=2.5m, with gown clips. Ergonomic, comfortable and foam bolstered headlight (grey and black), easy adjustable with two central adjustment wheels.



Modern focusable fiber optic headlight-system

with high performance optical system for true daylight colour rendering, wide focus range and powerful light spot.

Lightweight and comfortable headband, autoclavable adjustment stick. Bifurcated cable with gown clips for best use.

The Xenon Headlight System: composed of a 300W Xenon high power light source for advanced pure daylight, a fiber optic surgery headlight and a mobile trolley.

► **System I: Item code: 12-795C**

► **Item code: 08-595-S**

Fiber optic headlight

► **Item code: 12-204** End fitting Storz

► **Item code: 12-500X** Xenon 300W high power light source: Advanced Xenon lamp (5600°K) with sunlight spectrum for true color fidelity, light intensity adjustable by diaphragm, focus stepless adjustable, Power supply: 110/115V or 230/240V

► **Item code: 12-685**

Trolley

► **System II: Item code: 12-790C**

► **Item code: 08-590-S**

Fiber optic headlight

► **Item code: 12-204** End fitting Storz

► **Item code: 12-500X** Xenon 300W high power light source: Advanced Xenon lamp (5600°K) with sunlight spectrum for true color fidelity, light intensity adjustable by diaphragm, focus stepless adjustable, Power supply: 110/115V or 230/240V

► **Item code: 12-685**

Trolley

► **System III: Item code: 12-785C**

► **Item code: 08-585-S**

Fiber optic headlight

► **Item code: 12-204** End fitting Storz

► **Item code: 12-500X** Xenon 300W high power light source: Advanced Xenon lamp (5600°K) with sunlight spectrum for true color fidelity, light intensity adjustable by diaphragm, focus stepless adjustable, Power supply: 110/115V or 230/240V

► **Item code: 12-680**

Trolley

► Item code: **12-500X**

► Item code: **12-685**



► **System IV: Item code: 12-780C**

► Item code: **08-580-S**

Fiber optic headlight

► Item code: **12-204** End fitting Storz

► Item code: **12-500X** Xenon 300W

high power light source:

Advanced Xenon lamp (5600° K) with sunlight spectrum for true color fidelity, light intensity adjustable by diaphragm, focus stepless adjustable,

Power supply: 110/115V or 230/240V

► Item code: **12-680**

Trolley

Item code	Description	Headband
► 08-595-S	Headlight, super lightweight special cable	Black headband with top, two adjustment wheels
► 08-595-GR	Headlight, super lightweight special cable	Grey headband with top, two adjustment wheels
► 08-590-S	Headlight, super lightweight special cable	Black, open headband, one adjustment headband
► 08-590-GR	Headlight, super lightweight special cable	Grey, open headband, one adjustment headband
► 08-585-S	Headlight, shielded, light-weight special cable	Black headband with top, two adjustment wheels
► 08-585-GR	Headlight, shielded, light-weight special cable	Grey headband with top, two adjustment wheels
► 08-580-S	Headlight, shielded, light-weight special cable	Black, one headband, one adjustment wheel
► 08-580-GR	Headlight, shielded, light-weight special cable	Grey, one headband, one adjustment wheel

Light source end fittings (please order additionally):

► End fitting Storz Item code: **12-204**

► End fitting Wolf Item code: **12-202L**

Other end fittings on request.

Classification: Class I, according to appendix IX, RL 93/42 EWG

Spare parts:

► Item code: **99-X130P** Replacement bulb 300W

► Item code: **08-595K** Spare cable for Item code: **08-595** and **08-590**

► Item code: **08-585K** Spare cable for Item code: **08-585** and **08-580**

► Item code: **08-365R-BL** Headband (black)

► Item code: **08-365R-GR** Headband (grey)

► Item code: **08-365P** Spare neoprene

Spot size diameter

Working distance	20 cm	30 cm	60 cm
Spot size	ø 6-60 mm	ø 8-95 mm	ø 20-200 mm

► **FARO-CAM Headlights**

Micro Camera System: Item code: 12-0500V

Modern, lightweight high class solution for medical video documentation

Miniature camera

Videoprocessor, Videosystem PAL, NTSC

Videoresolution: 470 TV-lines

Shutter: max. 1/2.000.000 sec.

Video output 1xFBAS, 1xY/C(S-VHS)

Sensitivity 3 lux

Cable length 2,2m (longer on request)

► **Camera holder for Faromed headlights:**

Item code: 08-335MK-10



Video Documentation System

Modern High-End Solution for digital recording and archiving of video findings and diagnosis. Optional use with touch screen and foot switch „Made in Germany“.

► **Media Workstation**

Select: Item code: 19-100

Basic Edition: highly efficient functionality for recording, playing and archiving videos and individual images, database for an effective administration



Recording, Snapshot, Editing

► **Media Workstation**

Standard: Item code: 19-200

Comfortable and easy to use database for user, findings and diagnosis. Powerful tools for reports, editing videos and single picture, import of videos, picture and reports.

Languages: English, Spanish, German, French, Italian, Russian, Chinese (Mainland)



Digital Archiving

► **Foot-operated switch: Item code: 19-310**

► **Complete system on request**



Reports and Printouts